



张家港麦迪卡医疗器械有限公司
ZHANGJIAGANG MEDI MEDICAL EQUIPMENT CO., LTD.



CE ISO9001 ISO13485



Model: YA-YC HYDRAULIC OPERATION TABLE

Size: (2000±50mm)*(500±20mm) *(700-1000)±50mm
Application: For various surgical operations(head, neck etc.) also orthopedic operations with traction frame
Functions: 1.HI-LOW:(700-1000)±50mm 2.Head board fold upward≥30°and fold downward≥90° 3.Back board fold upward≥60°and fold downward≥20° 4.Leg board fold downward≥90° 5.Trendelenburg≥15°and Reverse-Trendelenburg≥15° 6.Left and Right lateral tilting angle≥15° 7.Head and leg section exchange to make complete examination in C-arm, X-ray photography
Structure: Steel structure
Accessories: 1 pair shoulder support, 1 pair armrests, 1 pair waist support frame , 1 pair leg support holder, 1 pc anesthesia screen frame, Hydraulic oil pump for elevation
Loading capacity: 180kgs
Packing: 133*76*97cm=0.98M3 G.W.:288kgs